

New England Biolabs Product Specification

Product Name: NEBNext[®] Ultra[™] II RNA Library Prep Kit for Illumina[®]
Catalog #: E7770S/L
Kit Components: NEBNext[®] Ultra[™] II End Prep Enzyme Mix (E7646)
NEBNext[®] First Strand Synthesis Enzyme Mix (E7761)
NEBNext[®] Second Strand Synthesis Enzyme Mix (E7425)
NEBNext[®] Ligation Enhancer (E7374)
NEBNext[®] Ultra[™] II End Prep Reaction Buffer (E7647)
NEBNext[®] First Strand Synthesis Reaction Buffer (E7421)
NEBNext[®] Second Strand Synthesis Reaction Buffer (E7427)
NEBNext[®] Ultra[™] II Ligation Master Mix (E7648)
NEBNext[®] Ultra[™] II Q5[®] Master Mix (E7649)
NEBNext[®] Adaptor Dilution Buffer (E7762)
0.1X TE (E7763)
Random Primers (E7422)
Nuclease Free Water (E7764)

Shelf Life: 18 months
Storage Temp: -20°C
Specification Version: PS-E7770S/L v2.0
Effective Date: 08 Apr 2022

Assay Name/Specification (minimum release criteria)

Functional Testing (Library Construction, RNA) - Each set of reagents is functionally validated and compared to the previous lot through construction of libraries made from commercially available RNA, using the kit's minimum and maximum input requirements. Libraries made from the previous and current lots for both input RNA amounts are sequenced together on the same Illumina flow cell and compared across various metrics including library yield, individual transcript abundance correlations (low vs. high input, old lot vs. new lot), 5'-3' transcript coverage, and fraction of reads mapping to a reference.

* **Individual Product Component Note** - Standard Quality Control Tests are performed for each component included in NEBNext[®] Ultra[™] II RNA Library Prep Kit for Illumina[®] and meet the designated specifications.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 08 Apr 2022

Derek Robinson
Director, Quality Control

