
New England Biolabs Product Specification

Product Name: *microRNA Marker*
Catalog #: *N2102S*
Shelf Life: *24 months*
Storage Temp: *-20°C*
Storage Conditions: *Provided in a ready-to-load solution containing 4M Urea and 0.04% Orange G*
Specification Version: *PS-N2102S v1.0*
Effective Date: *08 May 2019*

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Marker) - The banding pattern of microRNA Marker on a 15% denaturing TBE-Urea gel shows discrete, clearly identifiable bands at each band of the marker, when stained with SYBR Gold at a concentration of 1 X.



Date 08 May 2019

Derek Robinson
Director of Quality Control

