

New England Biolabs Product Specification

Product Name: Oligo d(T)₂₃ VN
Catalog #: S1327S
Shelf Life: 36 months
Storage Temp: -20°C
Storage Conditions: Supplied as dried, triethylammonium salt
Specification Version: PS-S1327S v2.0
Effective Date: 08 Nov 2022

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Oligonucleotide) - The banding pattern of Oligo d(T)₂₃ VN on a 15% TBE-Urea gel is evaluated against a control lot for expected size, sharpness, and relative intensity.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 08 Nov 2022

Lauren Brown
Quality Approver

