
New England Biolabs Product Specification

Product Name: Monarch[®] RNA Cleanup Columns (500 µg)
Catalog #: T2057L
Shelf Life: 18 months
Storage Temp: 25°C
Specification Version: PS-T2057L v2.0
Effective Date: 19 Jan 2021

Assay Name/Specification (minimum release criteria)

Functional Testing (RNA Recovery and Purity) - Twenty-four Monarch [®] RNA Cleanup Columns (T2057, 500 µg) are tested with 50 µg of RNA resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 19 Jan 2021

Derek Robinson
Director, Quality Control

