

pKLAC1

pKLAC1 is an expression vector capable both of replication in *E. coli* and stable integration into the genome of the yeast *Kluyveromyces lactis* (1). It is designed for high-level expression of recombinant protein in *K. lactis* using the *K. lactis* Protein Expression Kit (NEB #E1000).

In *E. coli*, it replicates using the pMB1 origin of replication from pBR322 (although the *rop* gene is missing) and carries the *bla* (Ap^R) marker for selection with ampicillin. Upon transformation of *K. lactis* GG799 competent cells (NEB #C1001), SacII- or BstXI-linearized pKLAC1 integrates into the *K. lactis* chromosome at the *LAC4* locus. Yeast transformants can be selected using the acetamidase selectable marker (*amdS*), which is expressed from the yeast *ADH1* promoter. Acetamidase expressed from pKLAC1 permits transformed cells to utilize acetamide as a sole nitrogen source on defined medium (2).

The multiple cloning site (MCS) is positioned to allow translational fusion of the *K. lactis* α -mating factor secretion domain (α -MF) to the N-terminus of the recombinant target protein. This directs the fusion protein to the general secretory pathway, but the α -MF domain is cleaved off in the Golgi apparatus by the Kex protease, resulting in secretion of the recombinant protein alone.

Expression of the recombinant fusion protein is driven by the *K. lactis* *LAC4* promoter, which has been modified to be transcriptionally silent in *E. coli* (1). This facilitates the cloning of proteins that are toxic to *E. coli*. This promoter is split such that when pKLAC1 is cleaved with SacII or BstXI, the recombinant protein and selectable marker are flanked by the two halves of the promoter. When these ends recombine with the *LAC4* promoter in the *K. lactis* chromosome, the result is integration of the recombinant fusion protein (driven by the *LAC4* promoter) and *amdS* upstream of the *LAC4* gene (driven by a duplicate copy of the *LAC4* promoter) (2).

Enzymes with unique restriction sites are shown in **bold** type and selected enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes can


be found on the NEB web site (choose Technical Reference > DNA Sequences and Maps). Restriction site coordinates refer to the position of the 5'-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Components of coordinated regions are indented below the region itself.

pMB1 origin of replication coordinates include the region from the –35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. Promoter and transcription terminator coordinates represent cloned regions and not necessarily the precise functional elements.

Feature	Coordinates	Source
expression region:		
α -mating factor		
leader sequence	14-334	<i>K. lactis</i>
MCS	257-343	–
<i>LAC4</i> TT region	355-937	<i>K. lactis</i>
<i>AdH1</i> promoter region	994-1696	<i>S. cerevisiae</i>
<i>amdS</i>	1697-3343	<i>A. nidulans</i>
<i>LAC4</i> promoter region (5' end)	4052-4632	<i>K. lactis</i>
origin	5674-5086	pMB1
<i>bla</i> (Ap ^R)	6705-5845	Tn3
<i>LAC4</i> promoter region (3' end)	7459-9091	<i>K. lactis</i> (modified)

ori = origin of replication
Ap = ampicillin
TT = transcription terminator



References

- (1) Colussi, P.A. and Taron, C.H. (2005) *Appl. Environ. Microbiol.*, 71, 7092–7098.
- (2) van Ooyen, A.J. et al. (2006) *FEMS Yeast Res.*, 6, 381–392.

GenBank Accession #: AY968582

See page 163 for ordering information.

There are no restriction sites for the following enzymes: AarI (x), AatII, AfeI, AflII, ApaI, AsiSI, BbvCI, BlnI, Bpu10I, BsiWI, FseI, FspAI (x), I-CeuI, I-SceI, MluI, MscI, PI-PspI, PI-SceI, PmeI, PmlI, PspOMI, PspXI, RsrII, SanDI (x), SfiI, SgrAI, SpeI, SrfI (x), SwaI, ZraI

(x) = enzyme not available from NEB